





Shenzhen Bioeasy Biotechnology Co., Ltd.

Add: A101, Building 2, Bioeasy Industry Park, No. 289 Yunchang Road,
Bao'an District, Shenzhen City, Guangdong Province, China
Tel: +86 755 2794 8546
Email: info@bioeasy.com
Web: www.bioeasy.com

Ver.202404



BIOEASY MILK SAFETY TEST

 Protect your milk from antibiotics, adulteration, mycotoxins and pesticides



01	Milk Antibiotics Screening	01
02	Milk Adulteration Screening	11
03	Aflatoxin M1 Screening in Milk	15
04	Milk Pesticides Screening	22
05	BIOEASY Instrument	24

Content

BIOEASY

BIOEASY founded in 2007, is a high-tech enterprise engaging in food safety, animal disease diagnosis and clinical diagnosis. Since the beginning, BIOEASY has constantly invested financial and material resources in new technology and innovation.

As an important part of milk safety, the contamination of antibiotics, adulteration, mycotoxins and pesticides in milk is a great concern to human health, which may cause allergy and hypersensitivity in human beings. They also pose a major problem to the dairy industry as antibiotic residues can inhibit or slow down starter culture in yogurt and cheese making, possibly resulting in major losses to the dairy concerned.

Working as a partner with farmers and processors, BIOEASY has kept pace with global testing requirements of the industry to ensure that milk and milk products are among the world safest foods!

Why BIOEASY

- 01 Total Solution**
A complete line of rapid tests, incubators and readers.
- 02 Meet Standards**
Meet EU, Russian, Australia and American standards.
- 03 Certificated**
ILVO, AOAC Research Institute Performance Tested MethodsSM
- 04 Professional R&D**
Positive and negative controls for quality insurance.



1

MILK ANTIBIOTICS SCREENING

Keeping Milk Free From Antibiotic Residues

Milk and other dairy products with drug residues beyond the MRLs, can pose serious health problems to consumers. Residues of antibiotics are mainly found in milk due to their injudicious usage in treating infectious diseases of animals. Moreover, some antibiotics are being used as feed additives indiscriminately which is another source of antibiotic residues in milk, ultimately responsible for potential public health importance.

Beta-lactams Rapid Test for Milk ▶


This rapid test is used for detection of Beta-lactams in milk based on the colloidal gold immunochromatography technology. It takes about 5 min for one test.




Advantages

- Rapid**
Result in 5 min
- Sensitive**
- User friendly**
Easily performed on site or in lab
- Visual interpretation**
- Cost-effective**
- Accurate and reliable**

Test Procedure



Results



Sensitivity

LOD: limits of detection (µg/kg – ppb)

Beta-lactams	BIOEASY LOD	EU MRL	Beta-lactams	BIOEASY LOD	EU MRL
Penicillin-G	1.5-2	4	Cephapirin	15-20	60
Ampicillin	3-4	4	Cefalonium	6-8	20
Amoxicillin	3-4	4	Cefazolin	40-50	50
Oxacillin	6-10	30	Cefoperazone	4-6	50
Cloxacillin	6-10	30	Cefquinome	12-18	20
Dicloxacillin	10-20	30	Ceftiofur	80-100	100
Nafcillin	20-30	30	Cefacetrile	25-30	125

Specification

Order code: YRM1001-2C Kit components: 96 tests/kit Shelf life: 18 months Storage: 2-8°C

*The information is for reference purpose only, please in kind prevail.


2IN1 BT(EU)
Beta-lactams + Tetracyclines
Rapid Test for Milk ▶

This rapid test is used for detection of Beta-lactams and Tetracyclines in milk based on the colloidal gold immunochromatography technology. It takes about 8 min for one test.


Advantages

- **Multiple**
2 in 1
- **Validated**
by ILVO
- **Rapid**
Result in 8 min
- **Sensitive**
- **User friendly**
Easily performed
on site or in lab
- **Visual interpretation**
- **Cost-effective**
- **Accurate and reliable**

Test Procedure



Results



Sensitivity

LOD: limits of detection (µg/kg – ppb)

Antibiotics	BIOEASY LOD	EU MRL	Antibiotics	BIOEASY LOD	EU MRL
Penicillin-G	1.5-2	4	Cephapirin	15-18	60
Ampicillin	3-4	4	Cefalonium	6-8	20
Amoxicillin	3-4	4	Cefazolin	40-50	50
Oxacillin	5-7	30	Cefoperazone	4-6	50
Cloxacillin	6-8	30	Cefquinome	12-18	20
Dicloxacillin	10-20	30	Ceftiofur	80-100	100
Nafcillin	20-30	30	Cefacetrile	25-30	125
Cefalexin	500-600	100	-	-	-
Doxycycline	30-50	n.f.u.	Oxytetracycline	30-50	100
Chlorotetracycline	30-50	100	Tetracycline	30-50	100

Specification

Order code:
YRM1007-40

Kit components:
96 tests/kit

Shelf life:
18 months

Storage:
2-8°C


2IN1 BTCef(EU)
Beta-lactams + Cefalexin + Tetracyclines
Rapid Test for Milk ▶

This rapid test is used for detection of Beta-lactams, Cefalexin and Tetracyclines in milk based on the colloidal gold immunochromatography technology. It takes about 9 min for one test.


Advantages

- **Multiple**
2 in 1
- **Validated**
by ILVO, AOAC
- **Rapid**
Result in 9 min
- **Sensitive**
- **User friendly**
Easily performed
on site or in lab
- **Visual interpretation**
- **Cost-effective**
- **Accurate and reliable**

Test Procedure



Results



Sensitivity

LOD: limits of detection (µg/kg – ppb)

Antibiotics	BIOEASY LOD	EU MRL	Antibiotics	BIOEASY LOD	EU MRL
Penicillin-G	1.5-2	4	Cefalexin	20-30	100
Ampicillin	3-4	4	Cephapirin	15-18	60
Amoxicillin	3-4	4	Cefalonium	6-8	20
Oxacillin	5-7	30	Cefazolin	40-50	50
Cloxacillin	6-8	30	Cefoperazone	4-6	50
Dicloxacillin	10-20	30	Cefquinome	12-18	20
Nafcillin	20-30	30	Ceftiofur	80-100	100
Cefacetrile	25-30	125	-	-	-
Doxycycline	30-50	n.f.u.	Oxytetracycline	30-50	100
Chlorotetracycline	30-50	100	Tetracycline	30-50	100

Specification

Order code:
YRM1008-40

Kit components:
96 tests/kit

Shelf life:
18 months

Storage:
2-8°C


3IN1 BST
Beta-lactams + Sulfonamides + Tetracyclines
Rapid Test for Milk ▶

This rapid test is used for detection of Beta-lactams, Sulfonamides, and Tetracyclines in milk based on the colloidal gold immunochromatography technology. It takes about 9 min for one test.


Advantages

- **Multiple**
3 in 1
- **Validated**
by ILVO, AOAC
- **Rapid**
Result in 9 min
- **Sensitive**
- **User friendly**
Easily performed
on site or in lab
- **Visual interpretation**
- **Cost-effective**
- **Accurate and reliable**

Test Procedure



Results



Sensitivity

LOD: limits of detection (µg/kg – ppb)

Antibiotics	LOD	Antibiotics	LOD	Antibiotics	LOD
Penicillin-G	1.5-2	Dicloxacillin	6-12	Cefazolin	50-70
Ampicillin	3-4	Piperacillin	6-12	Ceftiofur	80-100
Amoxicillin	3-4	Cefalexin	500-600	Cefquinome	12-18
Oxacillin	8-12	Cefacetrile	20-40	Cefoperazone	2-4
Nafcillin	20-35	Cephapirin	8-15	Cefalotin	30-50
Cloxacillin	6-10	Cefalonium	6-12	-	-
Sulfamethazine	15-25	Sulfadiazine	2-5	Sulfathiazole	2-5
Sulfapyridine	10-30	Sulfadimethoxypyrimidine	1.5-4	Sulfamethizole	2-5
Sulfaquinoxaline	14-28	Sulfadoxine	80-100	Sulfamonomethoxine	2-4
Sulfaphenazolum	20-40	Sulfamerazine	4-8	Sulfapyrazole	80-100
Sulfamethoxazole	3-6	Sulfamethoxypyridazine	20-40	Sulfamoxole	30-50
Sulfachloropyridazine	8-16	Sulfadimethoxine	12-24	Sulfisoxazole	80-100
Tetracycline	30-40	Oxytetracycline	30-40	-	-
Chlorotetracycline	30-40	Doxycycline	30-40	-	-

*The information is for reference purpose only, please in kind prevail.


AMINO 3IN1 GSN
Streptomycin + Gentamycin + Neomycin
Rapid Test for Milk ▶

This rapid test is used for detection of Streptomycin & Dihydrostreptomycin, Gentamycin, and Neomycin in milk based on the colloidal gold immunochromatography technology. It takes about 9 min for one test.


Advantages

- **Multiple**
3 in 1
- **Validated**
by ILVO
- **Rapid**
Result in 9 min
- **Sensitive**
- **User friendly**
Easily performed
on site or in lab
- **Visual interpretation**
- **Cost-effective**
- **Accurate and reliable**

Test Procedure



Results



Sensitivity

LOD: limits of detection (µg/kg – ppb)

Antibiotics	BIOEASY LOD	EU MRL	Antibiotics	BIOEASY LOD	EU MRL
Streptomycin	40-50	200	Gentamycin	30-40	100
Dihydrostreptomycin	40-50	200	Neomycin	80-100	1500

Specification

Order code: YRM1034
Kit components: 96 tests/kit
Shelf life: 18 months
Storage: 2-8°C

*The information is for reference purpose only, please in kind prevail.


4INI BSCT
Beta-lactams + Streptomycin + Chloramphenicol
+ Tetracyclines Rapid Test for Milk ▶

This rapid test is used for detection of Beta-lactams, Streptomycin, Chloramphenicol and Tetracyclines in milk based on the colloidal gold immunochromatography technology. It takes about 12 min for one test.


Advantages

- **Multiple**
4 in 1
- **Validated**
by ILVO
- **Rapid**
Result in 12 min
- **Sensitive**
- **User friendly**
Easily performed
on site or in lab
- **Visual interpretation**
- **Cost-effective**
- **Accurate and reliable**

Test Procedure



Results



Sensitivity

LOD: limits of detection (µg/kg – ppb)

Antibiotics	BIOEASY LOD	EU MRL	Antibiotics	BIOEASY LOD	EU MRL
Penicillin-G	1.5-2	4	Cefquinome	5-7	20
Ampicillin	2-3	4	Cefacetile	20-25	125
Amoxicillin	2-3	4	Cefalonium	3-5	20
Oxacillin	5-7	30	Cefazolin	15-20	50
Cloxacillin	3-6	30	Ceftiofur	80-100	100
Dicloxacillin	3-6	30	Cefoperazone	3-5	50
Nafcillin	8-10	30	Cephapirin	5-8	60
Doxycycline	7-10	n.f.u.	Oxytetracycline	7-10	100
Chlorotetracycline	7-10	100	Tetracycline	7-10	100
Streptomycin	80-100	200	Dihydrostreptomycin	50-60	200
Chloramphenicol	0.2-0.3	-	-	-	-

Specification

Order code:
YRM1039C

Kit components:
96 tests/kit

Shelf life:
18 months

Storage:
2-8°C


4INI BTSQ
Tetracyclines + Sulfonamides + Fluoroquinolones + Beta-lactams
Rapid Test for Milk ▶

This rapid test is used for detection of Tetracyclines, Sulfonamides, Fluoroquinolones and Beta-lactams in milk based on the colloidal gold immunochromatography technology. It takes about 9 min for one test.


Advantages

- **Multiple**
4 in 1
- **Validated**
by ILVO
- **Rapid**
Result in 9 min
- **Sensitive**
- **User friendly**
Easily performed
on site or in lab
- **Visual interpretation**
- **Cost-effective**
- **Accurate and reliable**

Test Procedure



Results



Sensitivity

LOD: limits of detection (µg/kg – ppb)

Antibiotics	LOD	Antibiotics	LOD	Antibiotics	LOD
Penicillin-G	1.5-2	Dicloxacillin	6-12	Cefazolin	50-70
Ampicillin	3-4	Piperacillin	6-12	Ceftiofur	80-100
Amoxicillin	3-4	Cefacetile	20-40	Cefquinome	12-18
Oxacillin	8-12	Cephapirin	8-15	Cefoperazone	2-4
Nafcillin	20-35	Cefalonium	6-12	Cefalotin	30-50
Cloxacillin	6-10	-	-	-	-
Sulfamethazine	15-25	Sulfadiazine	2-5	Sulfathiazole	2-5
Sulfapyridine	10-30	Sulfadimethoxypyrimidine	1.5-4	Sulfamethizole	2-5
Sulfaquinoxaline	14-28	Sulfadoxine	80-100	Sulfamonomethoxine	2-4
Sulfaphenazolum	20-40	Sulfamerazine	4-8	Sulfapyrazole	80-100
Sulfamethoxazole	3-6	Sulfamethoxypyridazine	20-40	Sulfamoxole	30-50
Sulfachloropyridazine	8-16	Sulfadimethoxine	12-24	Sulfisoxazole	80-100
Tetracycline	10-15	Oxytetracycline	5-10	-	-
Chlorotetracycline	20-30	Doxycycline	15-20	-	-
Enrofloxacin	10-15	Oxolinic Acid	50-70	Lomefloxacin	12-20
Ciprofloxacin	10-15	Cinoxacin	50-100	Sarafloxacin	7-10
Ofloxacin	6-10	Flumequine	50-70	Marbofloxacin	10-20
Pefloxacin	6-10	Danofloxacin	25-30	Orbifloxacin	40-50
Difloxacin	5-10	Norfloxacin	6-10	-	-

Specification

Order code:
YRM1059


Kit components:
96 tests/kit

Shelf life:
18 months

Storage:
2-8°C

16IN1
Rapid Test for Milk ▶

This 16in1 rapid test kit is based on the colloidal gold immunochromatography technology. It can be used for the detection of Monensin, Bacitracin, Spiramycin, Florfenicol, Trimethoprim, Colistin, Albendazole, Sulfonamides, Neomycin, Kanamycin, Spectinomycin, Gentamicin, Erythrocine, Lincomycin, Tylosin & Tilmicosin, Fluoroquinolones. It takes about 10 min for one test.



Advantages

Rapid
Result in 10 min


Multiple
16 in 1

User friendly
Easily performed
on site or in lab


Visual interpretation

Accurate and reliable

Test Procedure



Results



Sensitivity

	Channel 1 FMCS	Channel 2 SABT	Channel 3 NKSG	Channel 4 ELTF
T4	Spiramycin	Trimethoprim	Gentamicin	Fluoroquinolones
T3	Colistin	Bacitracin	Spectinomycin	Tylosin / Tilmicosin
T2	Monensin	Albendazole	Kanamycin	Lincomycin
T1	Florfenicol	Sulfonamides	Neomycin	Erythrocine

Specification

Order code:
YRM1112

Kit components:
16 tests/kit

Shelf life:
18 months

Storage:
2-8°C

BioQuick One step Rapid Test

BIOEASY BioQuick One Step Rapid Test is a new series of rapid test product, which can result by one step procedure and even in less time. With entirely new-looking, you can hold the cassette easier to add milk sample and interpret the result.



Advantages

Rapid
result in less time


One step
simplified procedure

User friendly
new outlook and easy
to hold

High specificity

Accurate and reliable

Test Procedure



Specification

Shelf life: 18 months **Storage:** 2-8°C

Art. No.	Products	Spec.	Test method
YR1M4001	BioQuick Test - 2IN1 BT Beta-lactams + Tetracyclines	30T & 100T	3 min 50±2°C
YR1M4002	BioQuick Test - 3IN1 BTCef Beta-lactams + Tetracyclines + Cefalexin	30T & 100T	3 min 50±2°C
YR1M4003	BioQuick Test - 4IN1 BTCC Beta-lactams + Tetracyclines + Ceftiofur + Cefalexin	30T & 100T	3 min 50±2°C
YR1M4004	BioQuick Test - 3IN1 BTCTf Beta-lactams + Tetracyclines + Ceftiofur	30T & 100T	3 min 50±2°C
YR1M4005	BioQuick Test - 4IN1 BTCC(RU) Beta-lactams + Tetracyclines + Ceftiofur + Cefalexin	30T & 100T	5 min 50±2°C
YR1M4006	BioQuick Test - 4IN1 BSTCTf Beta-lactams + Sulfonamides + Tetracyclines + Ceftiofur	30T & 100T	4 min 50±2°C
YR1M4007	BioQuick Test - 3IN1 BCC Beta-lactams + Cefalexin + Ceftiofur	30T & 100T	3 min 50±2°C
YR1M4008	BioQuick Test - AMINO 4IN1 GSKN Gentamycin + Streptomycin/Dihydrostreptomycin + Kanamycin + Neomycin	30T & 100T	4 min 40±2°C
YR1M4009	BioQuick Test - MACRO 4IN1 SLTE Spiramycin + Lincomycin + Tylosin/Tilmicosin + Erythromycin	30T & 100T	4 min 40±2°C
YR1M4010	BioQuick Test - 4IN1 BTSQ Sulfonamides + Fluoroquinolones + Beta-lactams + Tetracyclines	30T & 100T	4 min 50±2°C

*Contact BIOEASY for more details.
*The information is for reference purpose only, please in kind prevail.

*Contact BIOEASY for more details.
*The information is for reference purpose only, please in kind prevail.


2

MILK ADULTERATION SCREENING

Milk and dairy product adulteration came into global concern after breakthrough of melamine contamination. Adulterants in milk mainly include addition of vegetable protein, milk from different species, addition of whey and watering which are known as economically motivated adulteration. These adulterations do not pose any severe health risk. However, some adulterants (like melamine) are too harmful to be overlooked.

Rapid Test for Cow Milk Adulteration in Goat Milk and Sheep Milk

This rapid test is used for detection of cow milk in goat milk and sheep milk based on the colloidal gold immunochromatography technology. It takes about 6 min for one test.



Advantages



Rapid
Result in 6 min



Sensitive



User friendly
Easily performed on site or in lab



Visual interpretation




Cost-effective




Accurate and reliable

Test Procedure



Results



Specification

Shelf life: 18 months **Storage:** 2-8°C

Art. No.	Products	Sensitivity	Spec.	Test method
YREC1001-4-Y	Cow Milk Adulteration in Goat Milk and Sheep Milk	0.1-0.2%	96T	3+3 min 40±2°C
YREC1001-5-Y	Cow Milk Adulteration in Goat Milk and Sheep Milk	0.05-0.08%	96T	3+3 min 40±2°C

*The information is for reference purpose only, please in kind prevail.

Rapid Test for Cow Whey Adulteration in Milk

This rapid test is used for detection of cow whey composition in milk based on the colloidal gold immunochromatography technology. Without sample preparation required, it takes about 6 min for one test.

Advantages



Rapid
Result in 6 min



Sensitive



User friendly
Easily performed on site or in lab



Visual interpretation




Cost-effective




Accurate and reliable

Test Procedure



Results



Specification

Art. No.	Products	Sensitivity	Spec.	Test method
YREC1003	Cow Whey Adulteration in Milk	5%	48T	3+3 min 40±2°C
YREC1003-1	Cow Whey Adulteration in Milk	4-5%	48T	3+3 min 40±2°C

*The information is for reference purpose only, please in kind prevail.

Melamine Rapid Test for Milk

This rapid test is used for detection of Melamine in milk based on the colloidal gold immunochromatography technology.

Advantages



Rapid



Sensitive



User friendly
Easily performed on site or in lab



Visual interpretation




Cost-effective




Accurate and reliable

Test Procedure



Results



Specification

		Shelf life: 18 months Storage: 2-8°C		
Art. No.	Products	Sensitivity	Spec.	Test method
YRFM1003	Melamine Rapid Test for Dairy Products	100-500ppb	96T	6+7 min 40±2°C
YRM1003-10	Melamine Rapid Test for Milk	8-10ppb	96T	3+5 min 40±2°C
YRM1003-50	Melamine Rapid Test for Milk	30-50ppb	96T	3+3 min 40±2°C
YRM1003-100	Melamine Rapid Test for Milk	80-100ppb	96T	3+3 min 40±2°C
YRM1033-1	4IN1 MBCT (Melamine, Beta-lactams, Cephalaxin and Tetracyclines)	M: 30-50ppb	96T	3+9 min 40±2°C

*The information is for reference purpose only, please in kind prevail.

3



AFLATOXIN M1 SCREENING IN MILK

Milk and dairy products contamination with aflatoxin M1 is an important problem worldwide especially for developing countries for years. Because aflatoxins are carcinogenic to animals and perhaps humans, they are monitored closely in the food supply. The presence of aflatoxin M1 in these products is important issue, especially for children and infants, who are more susceptible than adults.



Aflatoxin M1 Rapid Test for Milk

This rapid test is used for detection of Aflatoxin M1 in milk and dairy products based on the colloidal gold immunochromatography technology.



Advantages



Rapid



Sensitive



User friendly
Easily performed on
site or in lab



Visual
interpretation



Cost-effective




Accurate and
reliable

Test Procedure




Add 200µL
milk and mix

Incubate 3 min
at 40±2°C



Insert
the test strip

Incubate 4 min
at 40±2°C




Visual interpretation



Interpretation by Reader

(Test procedure of YRM1004-3)

Results



Specification

Shelf life: 18 months Storage: 2-8°C

Art. No.	Products	Sensitivity	Spec.	Test method
YRFM1004-2	Aflatoxin M1 Rapid Test for Dairy Products	0.5ppb	96T	3+6 min 40±2°C
YRM1004-0.05	Aflatoxin M1 Rapid Test for Milk	0.04-0.05ppb	96T	2+8 min 40±2°C
YRM1004-0.5Q	Aflatoxin M1 Rapid Test for Milk	0.5ppb	96T	2+5 min 20-30°C
YRM1004-1	Aflatoxin M1 Rapid Test for Milk	0.08-0.1ppb	96T	3+4 min 40±2°C
YRM1004-3	Aflatoxin M1 Rapid Test for Milk	0.3-0.4ppb	96T	3+4 min 40±2°C
***	***	***	***	***

*The information is for reference purpose only, please in kind prevail.

Aflatoxin M1 Quantitative Rapid Test for Milk

This rapid test is used for quantitative detection of Aflatoxin M1 in milk based on the colloidal gold immunochromatography technology. It takes about 10 min for one test.

Advantages



Quantitative



Rapid
Result in 10 min



Sensitive



User friendly
Easily performed on site or in lab



Cost-effective



Accurate and reliable


Test Procedure



Add 200μL milk sample and mix well

Incubate 2 min at 40±2°C






Insert the test strip

Incubate 8 min at 40±2°C





Interpretation by Reader

(Test procedure of YRM1004D-2)

Tips: Please remember to insert the ID chip for each test kit in order to import the quantitative curve into the Reader.


Specification

		Shelf life: 18 months Storage: 2-8°C		
Art. No.	Products	Linear range	Spec.	Test method
YRM1004D-1	Aflatoxin M1 Quantitative Rapid Test for Milk	150-600 ppt	96T	2+8 min 40±2°C
YRM1004D-2	Aflatoxin M1 Quantitative Rapid Test for Milk	15-150 ppt	96T	2+8 min 40±2°C
YRM1004D-4	Aflatoxin M1 Quantitative Rapid Test for Milk	100-1000 ppt	96T	2+8 min 40±2°C


*The information is for reference purpose only, please in kind prevail.

How To Read Quantitative Result?


Tips: After the second incubation of the test strip, read the result directly within 1 min. The results are invalid after more than 1 min.




(Insert the ID chip)




01 Click "Test"




02 Select "Quantitative Test Strip"




03 Select the test item



04 Insert the test strip



05 Click "Test", automatic identification and detection



06 Automatically display the test result after a few seconds



Aflatoxin M1 Fluorescence Quantitative Rapid Test for Milk (15-150ppt)

This rapid test is used for quantitative detection of Aflatoxin M1 in milk based on fluorescence immunochromatography analysis technology. It takes about 12 min for one test. The linear range of this product is 15-150ppt.

Advantages

Quantitative

Rapid
Result in 12 min

Sensitive

User friendly
Easily performed on site or in lab

Cost-effective

Accurate and reliable

Test Procedure

Add 100μL of M1 sample diluent and 200μL of raw milk to microwells. Mix well

Incubate 2 min at 40 ± 2°C

Add 200μL of the mixture to sample wells

Incubate 10 min at 30±2°C

Interpretation by EASY-11 Analyzer

Tips: Please complete the reading of the result within 90 seconds after the end of the reaction.

Specification

		Shelf life: 18 months Storage: 2-30°C		
Art. No.	Product Name	Linear range	Spec.	Test method
YRM1004D-4K	Aflatoxin M1 Fluorescence Quantitative Rapid Test for Milk (15-150ppt)	15-150ppt	32T & 96T	2 min 40±2°C+10 min 30±2°C

*The information is for reference purpose only, please in kind prevail.

How to use EASY-11 Immunofluorescence Analyzer to read quantitative results?

Tips: Please remember to insert the ID chip for each new batch of test kit in order to upload the quantitative curve into the analyzer. Skip steps 2-3 if the batch curve already exists in the analyzer.

1 Log in

2 Insert the ID chip and click "Read ID Chip"

3 Load the curve successfully

4 Insert the test cassette into the test slot and click "Start Test"

5 Get the result (Please complete the reading of the result within 90 seconds after the end of the reaction)

Aflatoxin M1 ELISA Test Kit

This test kit is used for the detection of Aflatoxin M1 by the direct competitive ELISA method. The aflatoxin M1 in the sample competes with enzyme conjugate to bind the AFM1 antibody pre-coated in the microwell. TMB substrate is used for coloration and the absorbance value of the sample is negatively related to the content of aflatoxin M1, which is compared with the standard curve to obtain the corresponding residue content of Aflatoxin M1 in the sample.

Advantages



Quantitative



Sensitive



User friendly



Cost-effective




High quality



Accurate and reliable

Test Procedure



Specification

Shelf life: 12 months Storage: 2-8°C				
Art. No.	Products	Linear range	Spec.	Test method
YR201520	Aflatoxin M1 ELISA Test Kit	5-125ppt	96T	20+20+20 min 25°C

*The information is for reference purpose only, please in kind prevail.



MILK PESTICIDES SCREENING

In agriculture, crops like fruits, vegetables and cereals are treated with different types of pesticides. These chemicals accumulate in milk from contaminated grass, corn, silage and through direct application of pesticides to cattle. These pesticides cause a wide range of toxic effects and pose very severe health risks, specifically in infants with underdeveloped metabolic and enzymatic systems. Therefore, the screening of pesticide residues is very important to ensure the health and safety of consumers and realize national food safety.



5CE INSTRUMENT

Pesticide Residue Rapid Test for Milk

This rapid test is used for detection of pesticide residue in milk based on the colloidal gold immunochromatography technology.

Advantages



Rapid



Sensitive



User friendly
Easily performed on
site or in lab



Visual
interpretation




Cost-effective



Accurate and
reliable

Test Procedure



Specification

Shelf life: 18 months Storage: 2-8°C

Art. No.	Products	Sensitivity	Spec.	Test method
YRM1103	Carbaryl Rapid Test for Milk	20-25 ppb	96T	3+5 min 40±2°C
YRM1108	Imidacloprid Rapid Test for Milk	0.5-1 ppb	96T	3+3 min 40±2°C
...	...	...	...	...

*Contact BIOEASY for more details.

*The information is for reference purpose only, please in kind prevail.




YRM-IN03 Micro Heater



MiniT-N Incubator



YRM-IN05 Auto Incubator



YRM-IN04 BioQuick Incubator



EASY-2T Portable Incubator



YR-SG21 ReadEasy



OS-02 BioQuick Reader



EASY-11 Immunofluorescence Analyzer

*Contact BIOEASY for more details.

MILK SAFETY RAPID TEST						
No.	Art. No.	Products	Sensitivity(ppb)		Spec.	Test method
Adulteration						
1	YREC1001-4-Y	Rapid Test for Cow Milk Adulteration in Goat Milk and Sheep Milk	0.1-0.2%		96 T	3+3 min 40±2℃
2	YREC1001-5-Y	Rapid Test for Cow Milk Adulteration in Goat Milk and Sheep Milk	0.05-0.08%		96 T	3+3 min 40±2℃
3	YREC1003	Rapid Test for Cow Whey Adulteration in Milk	5%		48 T	3+3 min 40±2℃
4	YREC1003-1	Rapid Test for Cow Whey Adulteration in Milk	4-5%		48 T	3+3 min 40±2℃
5	YREC1003K	Fluorescence Semi-quantitative Rapid Test for Cow Whey Adulteration in Milk	10-15 ppm 15-30 ppm	30-50 ppm >50 ppm	30 T	8 min 30±2℃
Melamine						
6	YRFM1003	Melamine Rapid Test for Dairy Products	100-500		96 T	6+7 min 40±2℃
7	YRM1003-10	Melamine Rapid Test for Milk	8-10		96 T	3+5 min 40±2℃
8	YRM1003-50	Melamine Rapid Test for Milk	30-50		96 T	3+3 min 40±2℃
9	YRM1003-100	Melamine Rapid Test for Milk	80-100		96 T	3+3 min 40±2℃
Mycotoxin						
10	YRFM1004-2	Aflatoxin M1 Rapid Test for Dairy Products	0.5		96 T	3+6 min 40±2℃
11	YRFM1004-5	Aflatoxin M1 Rapid Test for Dairy Products	0.15-0.2	0.4-0.5	96 T	3+6 min 40±2℃
12	YRM1004-0.05	Aflatoxin M1 Rapid Test for Milk	0.04-0.05		96 T	2+8 min 40±2℃
13	YRM1004-0.05A	Aflatoxin M1 Rapid Test for Milk	0.04-0.05		96 T	2+8 min 40±2℃
14	YRM1004-0.5Q	Aflatoxin M1 Rapid Test for Milk	0.5		96 T	2+5 min 20-30℃
15	YRM1004-1	Aflatoxin M1 Rapid Test for Milk	0.08-0.1		96 T	3+4 min 40±2℃
16	YRM1004-3	Aflatoxin M1 Rapid Test for Milk	0.3-0.4		96 T	3+4 min 40±2℃
17	YRM1004D-1	Aflatoxin M1 Quantitative Rapid Test for Milk (150-600 ppt)	150-600 ppt		96 T	2+8 min 40±2℃
18	YRM1004D-2	Aflatoxin M1 Quantitative Rapid Test for Milk (15-150 ppt)	15-150 ppt		96 T	2+8 min 40±2℃
19	YRM1004D-4	Aflatoxin M1 Quantitative Rapid Test for Milk (100-1000 ppt)	100-1000 ppt		96 T	2+8 min 40±2℃
20	YRM1004D-4K	Aflatoxin M1 Fluorescence Quantitative Rapid Test for Milk (15-150 ppt)	15-150 ppt		32T & 96T	2 min 40±2℃+10 min 30±2℃
Hormone						
21	YRM1009	Dexamethasone & Prednisolone Rapid Test for Milk	Dexamethasone 0.3 Prednisolone 5-10		96 T	6+6 min 40±2℃
22	YRM1009-2	Dexamethasone & Betamethasone Rapid Test for Milk	Dexamethasone 0.15-0.2 Betamethasone 0.15-0.2 Prednisolone 0.4-0.8 Fludrocortisone 3-6 Hydrocortisone 15-30		96 T	3+5 min 40±2℃
23	YRM1053-1	Zeranol Rapid Test for Milk	0.8-1		96 T	6+4 min 40±2℃
24	YRM1140	Dexamethasone & Betamethasone Rapid Test for Milk	Dexamethasone 0.2-0.3 Betamethasone 0.2-0.3		96 T	3+5 min 40±2℃
25	YRM1142	Nandrolone Rapid Test for Milk	1-2		96 T	3+3 min 40±2℃
β-lactamase						
26	YRM1002	Beta-lactamase Rapid Test for Milk	0.35-0.4 IU/L 0.18-0.2 IU/L 0.12-0.14 IU/L		96 T	x+3+3 min 40±2℃
27	YRM1002-3	Beta-lactamase Rapid Test for Milk	3-4 U/mL		96 T	6+3+3 min 40±2℃
Antibiotic residue						
28	YRM1001-2C	Rapid Test for Beta-lactams in Milk	PG: 1.5-2...		96 T	2+3 min 40±2℃
29	YRM1001-3	Beta-lactams+Cefalexin Rapid Test for Milk	PG: 1.5-2...		96 T	3+5 min 40±2℃
30	YRM1005-10	Streptomycin Rapid Test for Milk	Streptomycin 10-20 Dihydrostreptomycin 10-20		96 T	3+3 min 40±2℃
31	YRM1006-20	Fluoroquinolones Rapid Test for Milk	Enrofloxacin: 3-5...		96 T	3+3 min 40±2℃

MILK SAFETY RAPID TEST					
No.	Art. No.	Products	Sensitivity(ppb)	Spec.	Test method
32	YRM1007-5	Gentamicin Rapid Test for Milk	5-6	96 T	3+3 min 40±2℃
33	YRM1007-40	2IN1 BT(EU) Beta-lactams+Tetracyclines Rapid Test for Milk	PG: 1.5-2...	96 T	3+5 min 40±2℃
34	YRM1007-40T	2IN1 BT(EU) Beta-lactams+Tetracyclines Rapid Test for Milk	PG: 2-3...	96 T	3+3 min 40±2℃
35	YRM1007-100	Gentamicin Rapid Test for Milk	60-100	96 T	3+5 min 40±2℃
36	YRM1008-1C	2IN1 BT(RU) Beta-lactams+Tetracyclines Rapid Test for Milk	PG: 1.5-2...	96 T	3+5 min 40±2℃
37	YRM1008-40	2IN1 BTCef(EU) Beta-lactams+Cefalexin+Tetracyclines Rapid Test for Milk	PG: 1.5-2...	96 T	3+6 min 40±2℃
38	YRM1010-8	Chloramphenicol Rapid Test for Milk	0.08-0.1	96 T	3+5 min 40±2℃
39	YRM1011-5	Tetracyclines Rapid Test for Milk	Chlortetracycline 4-5 Oxytetracycline 2-3 Doxycycline 4-5 Tetracycline 4-5	96 T	3+6 min 40±2℃
40	YRM1013-5	Lincomycin Rapid Test for Milk	3-4	96 T	1+5 min 40±2℃
41	YRM1013-100	Lincomycin Rapid Test for Milk	75-100	96 T	3+5 min 40±2℃
42	YRM1014-5	Kanamycin Rapid Test for Milk	5-10	96 T	3+3 min 40±2℃
43	YRM1014-50	Kanamycin Rapid Test for Milk	40-80	96 T	3+5 min 40±2℃
44	YRM1015	MACRO 2IN1 Tylosin&Tilmicosin Rapid Test for Milk	Tylosin 10-20 Tilmicosin 25-30	96 T	3+3 min 40±2℃
45	YRM1016A	2IN1 QS Quinolones and Sulfonamides Rapid Test for Milk	Enrofloxacin:3-5...	96 T	3+5 min 40±2℃
46	YRM1016C	2IN1 SQ Sulfonamides and Quinolones Rapid Test for Milk	Enrofloxacin: 1-2...	96 T	5+5 min 40±2℃
47	YRM1021	Aminoglycosides 3IN1 SKG (Streptomycin+Gentamycin+Kanamycin) Rapid Test Kit	Streptomycin 30-40 Dihydrostreptomycin 30-40 Gentamycin 5-10 Kanamycin 5-10	96 T	3+6 min 40±2℃
48	YRM1022	MACRO 3IN1 ELT Erythromycin+Lincomycin+Tylosin&Tilmicosin Rapid Test for Milk	Erythromycin 3-5 Lincomycin 3-6 Tylosin 10-20 Tilmicosin 40-50	96 T	3+7 min 40±2℃
49	YRM1024-3	Sulfonamides Rapid Test for Milk	Sulfamethazine: 8-12...	96 T	3+3 min 40±2℃
50	YRM1026	3IN1 BST Beta-lactams+Sulfonamides+Tetracyclines Rapid Test for Milk	PG: 1.5-2...	96 T	3+6 min 40±2℃
51	YRM1026-3	3IN1 Rapid Test for Milk (Beta-lactams+Tetracyclines+Sulfonamides+Ceftiofur)	PG: 1-2...	96 T	1+4 min 40±2℃
52	YRM1031-2	Trimethoprim Rapid Test for Milk	25-35	96 T	3+5 min 40±2℃
53	YRM1034	Amino 3IN1 GSN Streptomycin+Gentamycin+Neomycin Rapid Test for Milk	Streptomycin 40-50 Dihydrostreptomycin 40-50 Gentamycin 30-40 Neomycin 80-100	96 T	3+6 min 40±2℃
54	YRM1039C	4IN1 BSCT Beta-lactams+Streptomycin+Chloramphenicol+Tetracyclines Rapid Test for Milk	PG: 1.5-2...	96 T	5+7 min 40±2℃
55	YRM1039C-10	4IN1 BSCT Beta-lactams+Streptomycin+Chloramphenicol+Tetracyclines Rapid Test for Milk	PG: 1.5-2...	96 T & 48T	3+7 min 40±2℃
56	YRM1039R	4IN1 BSCT Beta-lactams+Streptomycin+Chloramphenicol+Tetracyclines Rapid Test for Milk	PG: 2-3...	96 T	3+7 min 40±2℃
57	YRM1039CQ	4IN1 BCTC Beta-lactams+Ceftiofur+Tetracyclines+Chloramphenicol Rapid Test for Milk	PG: 2-3...	96 T	3+7 min 40±2℃
58	YRM1054-1	Florfenicol & Thiamphenicol Rapid Test for Milk	Florfenicol 0.08-0.1 Thiamphenicol 0.5-0.8	96 T	3+3 min 40±2℃
59	YRM1059	4IN1 BTSQ Tetracyclines+Sulfonamides+Fluoroquinolones+Beta-lactams Rapid Test for Milk	Tetracycline: 10-15...	96 T	3+6min 40±2℃
60	YRM1060	Colistin Rapid Test for Milk	8-10	96 T	3+3min 40±2℃
61	YRM1062	Rapid Test for Bacitracin in Milk	15-20	96 T	5+5min 40±2℃
62	YRM1065	Monensin Rapid Test for Milk	25-35	96 T	3+3min 40±2℃
63	YRM1065-1	Monensin Rapid Test for Milk	3-5	96 T	3+5min 40±2℃
64	YRM1074	MACRO 3IN1(ETS) Erythromycin+Tylosin/Tilmicosin+Spiramycin Rapid Test for Milk	Erythromycin 15-20 Tylosin 25-30 Spiramycin 80-100 Tilmicosin 25-30	96 T	3+6min 40±2℃

*The information is for reference purpose only, please in kind prevail.

MILK SAFETY RAPID TEST					
No.	Art. No.	Products	Sensitivity(ppb)	Spec.	Test method
65	YRM1075	3IN1(QLT) Quinolones+Lincomycin+ Tylosin Rapid Test for Milk	Tylosin: 20-30...	96 T	3+5min 40±2°C
66	YRM1076	Apramycin Rapid Test for Milk	15-20	96 T	3+3min 40±2°C
67	YRM1077	Spiramycin Rapid Test for Milk	15-20	96 T	3+5min 40±2°C
68	YRM1080	4IN1(RU) NKSG Neomycin+Kanamycin+ Spectinomycin + Gentamicin Rapid Test for Milk	Neomycin 40-60 Kanamycin 40-50 Spectinomycin 5-10 Gentamicin 10-20	96 T	3+7min 40±2°C
69	YRM1081	4IN1 ELTF Erythromycin + Lincomycin + Tylosin/Timicosin + Fluoroquinolones Rapid Test for Milk	Erythromycin: 4-6...	96 T	3+7min 40±2°C
70	YRM1092	5IN1 CBSCT Cefalexin+Beta-lactams+Streptomycin+ Chloramphenicol+Tetracyclines Rapid Test for Milk	PG: 2-3...	96 T	3+7min 40±2°C
71	YRM1098	4IN1 CGSN Chloramphenicol + Gentamycin + Streptomycin +Neomycin Rapid Test for Milk	Chloramphenicol 0.2-0.3 Gentamycin 15-20 Streptomycin 30-50 Dihydrostreptomycin 30-50 Neomycin 60-80	96 T	3+6 min 40±2°C
72	YRM1105	Pirlimycin Rapid Test for Milk	40-50	96 T	3+5 min 40±2°C
73	YRM1105-1	Pirlimycin Rapid Test for Milk	150-200	96 T	3+5 min 40±2°C
74	YRM1112	16in1 Rapid Test for Milk	Florfenicol: 0.2-0.3	16 T	3+7 min 40±2°C
75	YRM1117	Rifaximin Rapid Test for Milk	Rifaximin 4-6 Rifampin 8-15	96 T	3+5 min 40±2°C
76	YRM1118	Diclofenac Rapid Test for Milk	3-5	96 T	3+6 min 40±2°C
77	YRM1121	AMINO 5IN1 SKGSN Spectinomycin + Kanamycin + Gentamycin + Streptomycin/Dihydrostreptomycin + Neomycin Rapid Test for Milk	Spectinomycin 4-8 Kanamycin 30-60 Gentamycin 4-8 Streptomycin 15-30 Dihydrostreptomycin 15-30 Neomycin 10-20	96 T	3+5 min 40±2°C
78	YRM1125	2IN1 TS Trimethoprim+Sulfonamides Rapid Test for Milk	Trimethoprim: 10-15...	96 T	3+5 min 40±2°C
79	YRM1126	4IN1 TCCB Tetracyclines+ Ceftiofur+ Cephalexin+ Beta-lactams Rapid Test for Milk	PG: 0.5-1.0...	96 T	1+5 min 40±2°C
80	YRM1127	2IN1 BCTf Beta-lactams + Ceftiofur Rapid Test for Milk	PG: 1-2...	96 T	1+4 min 40±2°C
81	YRM1128	4IN1 RPNF Rifaximin+ Pirlimycin+ Novobiocin+ Flavomycin Rapid Test for Milk	Rifaximin 5-10 Rifampicin 10-30 Pirlimycin 50-100 Clindamycin 2-5 Lincomycin 10-20 Novobiocin 20-50 Flavomycin 150-200	96 T	3+7 min 40±2°C
82	YRM1129	AMINO 2IN1 Neomycin+ Spectinomycin Rapid Test for Milk	Neomycin 30-60 Spectinomycin 20-30	96 T	3+5 min 40±2°C
83	YRM1130	MACRO 4IN1 SLTE Spiramycin + Lincomycin + Tylosin/Tilmicosin + Erythromycin Rapid Test for Milk	Spiramycin 4-8 Lincomycin 6-12 Tylosin 2-4 Tilmicosin 4-5 Erythromycin 3-4	96 T	3+5 min 40±2°C
84	YRM1131	3IN1 BCC(EU) Beta-lactams+Cefalexin+ Ceftiofur Rapid Test for Milk	PG: 1.5-2...	96 T	3+6 min 40±2°C
85	YRM1133	5IN1 CBCTS Ceftiofur+ Beta-lactams+ Cefalexin+ Tetracyclines+ Sulfonamides Rapid Test for Milk	PG: 1-2...	96 T	3+7 min 40±2°C
86	YRM1138	Abamectins Rapid Test for Milk	Abamectin 0.8-1.0 Ivermectin 1.0-1.5 Eprinomectin 1.5-3.0 Doramectin 12-14	48 T	3+6 min 10-30°C
87	YRM1146	4IN1 BSTC Beta-lactams+Sulfonamides+Tetracyclines+ Chloramphenicol Rapid Test for Milk	PG: 1.5-2...	96 T	3+6 min 40±2°C
88	YRM1147	5IN1 LETFS Lincomycin+Erythromycin+Tylosin/Tilmicosin+ Fluoroquinolones+Spiramycin Rapid Test for Milk	Lincomycin: 3-5...	96 T	3+7 min 40±2°C
89	YRM1148	5IN1 NKSGS Neomycin+Kanamycin+Spectinomycin+ Gentamicin+Sulfonamides Rapid Test for Milk	Neomycin: 40-60...	96 T	3+7 min 40±2°C
90	YRM1149	4IN1 GLSE Gentamicin+Lincomycin+Streptomycin/ Dihydrostreptomycin+Erythromycin Rapid Test for Milk	Gentamicin 5-10 Lincomycin 6-12 Streptomycin 30-50 Dihydrostreptomycin 30-50 Erythromycin 3-6	96 T	3+5 min 40±2°C

MILK SAFETY RAPID TEST					
No.	Art. No.	Products	Sensitivity(ppb)	Spec.	Test method
91	YRM1154M	2IN1 FC Florfenicol&Thiamphenicol+Chloramphenicol Rapid Test for Milk	Florfenicol 0.1 Thiamphenicol 0.8 Chloramphenicol 0.1	96 T	10+5 min 40±2°C
Veterinary drug residue					
92	YRM1052-1	Metamizole metabolite(4-Aminoantipyrine) Rapid Test for Milk	10-20	96 T	3+3 min 40±2°C
93	YRM1052M	Sulpyrine/Dipyrone/Metamizole Rapid Test for Milk	Metamizole and its metabolite 30-50	96 T	3+3 min 40±2°C
94	YRM1061	Diminazene Rapid Test for Milk	15-20	96 T	3+3 min 40±2°C
95	YRM1063	Flunixin Rapid Test for Milk	25-35	96 T	3+3 min 40±2°C
96	YRM1063-1	Flunixin Rapid Test for Milk	5-Hydroxyflunixin 20-30 Flunixin 2-4	96 T	3+3 min 40±2°C
97	YRM1063-2	Flunixin and its Metabolites Rapid Test for Milk	5-Hydroxyflunixin 15-20 Flunixin 2-4	96 T	3+3 min 40±2°C
98	YRM1064	Albendazole Rapid Test for Milk	30-50	96 T	3+5 min 40±2°C
99	YRM1066	Clopidol Rapid Test for Milk	10-14	48 T	3+6 min 40±2°C
100	YRM1066-3	Clopidol Rapid Test for Milk	2-3	96 T	3+6 min 40±2°C
101	YRM1115	Meloxicam Rapid Test for Milk	1-2	96 T	3+6 min 40±2°C
102	YRM1120	Halofuginone Rapid Test for Milk	6-10	96 T	3+3 min 40±2°C
103	YRM1122	Levamisole Rapid Test for Milk	3-5	96 T	3+4 min 40±2°C
104	YRM1124	Isometamidium Rapid Test for Milk	40-60	96 T	3+3 min 40±2°C
105	YRM1143	Toltrazuril and its metabolites Rapid Test for Milk	Toltrazuril 6-10 Toltrazuril sulfone 3-6 Toltrazuril sulfoxide 2-4	96 T	5+3 min 40±2°C
106	YRM1144	Carprofen Rapid Test for Milk	5-10	96 T	3+3 min 40±2°C
107	YRM1145	2IN1 MF Meloxicam+Flunixin Rapid Test for Milk	Meloxicam 1.5-3 Flunixin 4-8 5-Hydroxyflunixin 15-20	96 T	5+5 min 40±2°C
Pesticide residue					
108	YRM1103	Carbaryl Rapid Test for Milk	20-25	96 T	3+5 min 40±2°C
109	YRM1108	Imidacloprid Rapid Test for Milk	0.5-1	96 T	3+3 min 40±2°C
110	YRM1109	Paraquat Rapid Test for Milk	1-2	96 T	3+3 min 40±2°C
111	YRM1110	Acetamiprid Rapid Test for Milk	1-1.5	96 T	3+3 min 40±2°C
112	YRM1113	Chlorothalonil Rapid Test for Milk	10-15	96 T	3+5 min 40±2°C
113	YRM1114	Cyromazine Rapid Test for Milk	5-8	96 T	3+3 min 40±2°C
114	YRM1116	Fipronil Rapid Test for Milk	15-20	96 T	3+5 min 40±2°C
115	YRM1119	Thiamethoxam Rapid Test for Milk	5-10	96 T	3+5 min 40±2°C
116	YRM1136	Sulfoxaflor Rapid Test for Milk	5-10	96 T	3+5 min 40±2°C
117	YRM1137	Carbofuran Rapid Test for Milk	2-3	96 T	3+5 min 40±2°C
		Others			
118	YRM1033-1	4IN1 MBCT Melamine+Beta-lactams+Cephalexin+ Tetracyclines Rapid Test for Milk	PG: 0.8-1.5...	96 T	3+9 min 40±2°C
119	YRM1102-1	Carbendazim & Albendazole Rapid Test for Milk	Carbendazim 20-25 Albendazole sulfone 8-12 Albendazole sulfoxide 8-12 Albendazole 30-50	96 T	3+5 min 40±2°C
120	YRM1102-4	Carbendazim & Albendazole Rapid Test for Milk	Carbendazim 2-4 Albendazole sulfone 4-8 Albendazole sulfoxide 4-8 Albendazole 12-14	96 T	3+5 min 40±2°C
121	YRM1139	Diazepam Rapid Test for Milk	4+6	96 T	3+6 min 40±2°C

*See page 10 for a catalog of BioQuick One-step Test Cassette.

*The information is for reference purpose only, please in kind prevail.